



#accessible

Conservative exploration due to micro-tip design
Anatomical 3-D orientation
Precision control

AXS_tome



Cold polyp resection
in a new dimension

Perfect resection thanks to unique tube design
Multiple snare diameters for diminutive polyps

**CrossSnare
ZERO**



Mono pull wire prevents the branches from tilting
Ensures large sizes of sample material

EndoBite²



#safe

Intuitive operation through All-in-One System
twist 'n' CATCH basket for fast and reliable retrieval

ROCKSTAR



#extreme

Precise and kink-resistant needle activation
40% lower injection resistance
Easy handling due to single hand operation

MANTA



emed
confidence in performance

Gastroenterology in mind
Precision in Endoscopic Procedures
Smart Device System
Automatic adjustment
Ergonomics and comfort
Built in Argon mode
Wireless Footswitch

